

Product datasheet for **GA112804**

Human TRPM3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRPM3 CRISPRa kit - CRISPR gene activation of human transient receptor potential cation channel subfamily M member 3
Symbol:	TRPM3
UniProt ID:	Q9HCF6
Synonyms:	GON-2; LTRPC3; MLSN2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)