

Product datasheet for **GA112807**

Human **PGAP1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	PGAP1 CRISPRa kit - CRISPR gene activation of human post-GPI attachment to proteins 1
Symbol:	PGAP1
UniProt ID:	Q75T13
Synonyms:	Bst1; ISPD3024; MRT42; SPG67
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)