

Product datasheet for **GA112821**

Human PIF1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIF1 CRISPRa kit - CRISPR gene activation of human PIF1 5'-to-3' DNA helicase
Symbol:	PIF1
Components:	GA112821G1, PIF1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112821G2, PIF1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112821G3, PIF1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H611</u>
Synonyms:	C15orf20; PIF
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)