

## Product datasheet for **GA112843**

### Human CWH43 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	CWH43 CRISPRa kit - CRISPR gene activation of human cell wall biogenesis 43 C-terminal homolog
Symbol:	CWH43
Components:	GA112843G1, CWH43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112843G2, CWH43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112843G3, CWH43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<a href="#">Q9H720</a>
Synonyms:	CWH43-C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)