

Product datasheet for **GA112867**

Human OPA3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	OPA3 CRISPRa kit - CRISPR gene activation of human outer mitochondrial membrane lipid metabolism regulator OPA3
Symbol:	OPA3
UniProt ID:	Q9H6K4
Synonyms:	MGA3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)