

Product datasheet for **GA112869**

Human PROSER1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PROSER1 CRISPRa kit - CRISPR gene activation of human proline and serine rich 1
Symbol:	PROSER1
Components:	GA112869G1, PROSER1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112869G2, PROSER1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112869G3, PROSER1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86XN7</u>
Synonyms:	C13orf23
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)