

Product datasheet for **GA112874**

Human **CFAP43** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CFAP43 CRISPRa kit - CRISPR gene activation of human cilia and flagella associated protein 43
Symbol:	CFAP43
Components:	GA112874G1, CFAP43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112874G2, CFAP43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112874G3, CFAP43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NDM7</u>
Synonyms:	bA373N18.2; C10orf79; WDR96
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

