

Product datasheet for **GA112881**

Human PAAF1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PAAF1 CRISPRa kit - CRISPR gene activation of human proteasomal ATPase associated factor 1
Symbol:	PAAF1
Components:	GA112881G1, PAAF1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112881G2, PAAF1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112881G3, PAAF1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BRP4</u>
Synonyms:	PAAF; Rpn14; WDR71
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)