


Product datasheet for **GA112886**

Human FP100 (FAAP100) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FAAP100 CRISPRa kit - CRISPR gene activation of human FA core complex associated protein 100
Symbol:	FAAP100
Components:	GA112886G1, FP100 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112886G2, FP100 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112886G3, FP100 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q0VG06</u>
Synonyms:	chromosome 17 open reading frame 70; FAAP100; Fanconi anemia associated protein 100 kDa subunit; Fanconi anemia core complex 100 kDa subunit; FLJ22175; FLJ30151
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)