

## Product datasheet for **GA112891**

### Human SLC35F5 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SLC35F5 CRISPRa kit - CRISPR gene activation of human solute carrier family 35 member F5
Symbol:	SLC35F5
Components:	GA112891G1, SLC35F5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112891G2, SLC35F5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112891G3, SLC35F5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<a href="#">Q8WV83</a>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)