

Product datasheet for **GA112894**

Human PHAF1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C16orf70 CRISPRa kit – CRISPR gene activation of human chromosome 16 open reading frame 70
Symbol:	PHAF1
Components:	GA112894G1, PHAF1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112894G2, PHAF1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112894G3, PHAF1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BSU1</u>
Synonyms:	C16orf6; lin-10; LIN10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

