

Product datasheet for **GA112899**

Human IP3KC (ITPKC) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ITPKC CRISPRa kit - CRISPR gene activation of human inositol-trisphosphate 3-kinase C
Symbol:	ITPKC
Components:	GA112899G1, IP3KC gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112899G2, IP3KC gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112899G3, IP3KC gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96DU7</u>
Synonyms:	IP3-3KC; IP3KC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)