

## Product datasheet for **GA112906**

### Human C2orf44 (WDCP) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	WDCP CRISPRa kit - CRISPR gene activation of human WD repeat and coiled coil containing
Symbol:	WDCP
Components:	GA112906G1, C2orf44 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112906G2, C2orf44 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112906G3, C2orf44 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q9H6R7</a></u>
Synonyms:	C2orf44; PP384
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)