

Product datasheet for **GA112914**

Human **LRRC27** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC27 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 27
Symbol:	LRRC27
Components:	GA112914G1, LRRC27 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112914G2, LRRC27 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112914G3, LRRC27 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9C0I9</u>
Synonyms:	leucine rich repeat containing 27; OTTHUMP00000020759; OTTHUMP00000020760
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)