

## Product datasheet for **GA112918**

### Human ZKSCAN3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	ZKSCAN3 CRISPRa kit - CRISPR gene activation of human zinc finger with KRAB and SCAN domains 3
Symbol:	ZKSCAN3
UniProt ID:	<u><a href="#">Q9BRR0</a></u>
Synonyms:	dj874C20.1; dj874C20.1.; ZF47; zfp-47; Zfp47; ZFP306; ZNF306; ZNF309; ZSCAN13; ZSCAN35
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)