

Product datasheet for **GA112931**

Human KCNIP4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNIP4 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel interacting protein 4
Symbol:	KCNIP4
UniProt ID:	Q6PIL6
Synonyms:	CALP; KCHIP4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)