


Product datasheet for **GA112937**

Human SEL1L2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SEL1L2 CRISPRa kit - CRISPR gene activation of human SEL1L2 adaptor subunit of ERAD E3 ligase
Symbol:	SEL1L2
Components:	GA112937G1, SEL1L2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112937G2, SEL1L2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112937G3, SEL1L2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5TEA6</u>
Synonyms:	C20orf50
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)