

Product datasheet for **GA112964**

Human **LY6G5C** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LY6G5C CRISPRa kit - CRISPR gene activation of human lymphocyte antigen 6 family member G5C
Symbol:	LY6G5C
Components:	GA112964G1, LY6G5C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112964G2, LY6G5C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112964G3, LY6G5C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5SRR4</u>
Synonyms:	C6orf20; G5c; G5c, NG33, C6orf20, LY6G5CA, LY6G5CB; LY6G5CA; LY6G5CB; lymphocyte antigen-6 G5C; lymphocyte antigen 6 complex, locus G5C; lymphocyte antigen 6 complex G5C; NG33
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)