

Product datasheet for **GA112988**

Human ZNF436 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF436 CRISPRa kit - CRISPR gene activation of human zinc finger protein 436
Symbol:	ZNF436
Components:	GA112988G1, ZNF436 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112988G2, ZNF436 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112988G3, ZNF436 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9C0F3
Synonyms:	Zfp46; ZNF
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)