

Product datasheet for **GA113000**

Human SH3BP5L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SH3BP5L CRISPRa kit - CRISPR gene activation of human SH3 binding domain protein 5 like
Symbol:	SH3BP5L
Components:	GA113000G1, SH3BP5L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113000G2, SH3BP5L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113000G3, SH3BP5L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q7L8J4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)