


Product datasheet for **GA113004**

Human LNP (LNPK) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LNPK CRISPRa kit - CRISPR gene activation of human lunapark, ER junction formation factor
Symbol:	LNPK
Components:	GA113004G1, LNP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113004G2, LNP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113004G3, LNP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9C0E8
Synonyms:	LNPK; LNP1; UI; ulnaless
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)