

Product datasheet for **GA113045**

Human SYNC activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SYNC CRISPRa kit - CRISPR gene activation of human syncoilin, intermediate filament protein
Symbol:	SYNC
UniProt ID:	Q9H7C4
Synonyms:	SYNC1; SYNCOILIN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)