

Product datasheet for **GA113057**

Human TDRD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TDRD3 CRISPRa kit - CRISPR gene activation of human tudor domain containing 3
Symbol:	TDRD3
UniProt ID:	Q9H7E2
Synonyms:	RP11-459E2.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)