

## Product datasheet for **GA113091**

### Human TMEM163 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TMEM163 CRISPRa kit - CRISPR gene activation of human transmembrane protein 163
Symbol:	TMEM163
UniProt ID:	<a href="#">Q8TC26</a>
Synonyms:	DC29; SV31
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)