

Product datasheet for **GA113093**

Human CAB39L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CAB39L CRISPRa kit - CRISPR gene activation of human calcium binding protein 39 like
Symbol:	CAB39L
UniProt ID:	<u>Q9H9S4</u>
Synonyms:	antigen MLAA-34; bA103J18.3; calcium binding protein 39-like; FLJ12577; MO2L; MO2L, FLJ12577, bA103J18.3, RP11-103J18.3; MO25-BETA; OTTHUMP00000018410; RP11-103J18.3; sarcoma antigen NY-SAR-79; U937-associated antigen
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)