

Product datasheet for **GA113094**

Human ITM2C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ITM2C CRISPRa kit - CRISPR gene activation of human integral membrane protein 2C
Symbol:	ITM2C
Components:	GA113094G1, ITM2C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113094G2, ITM2C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113094G3, ITM2C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NQX7</u>
Synonyms:	BRI3; BRICD2C; E25; E25C; ITM3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)