

Product datasheet for **GA113095**

Human TSPAN14 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSPAN14 CRISPRa kit - CRISPR gene activation of human tetraspanin 14
Symbol:	TSPAN14
UniProt ID:	<u>Q8NG11</u>
Synonyms:	DC-TM4F2; TM4SF14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)