

Product datasheet for **GA113102**

Human C1orf25 (TRMT1L) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRMT1L CRISPRa kit - CRISPR gene activation of human tRNA methyltransferase 1 like
Symbol:	TRMT1L
Components:	GA113102G1, C1orf25 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113102G2, C1orf25 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113102G3, C1orf25 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z2T5</u>
Synonyms:	bG120K12.3; C1orf25; MST070; MSTP070; TRM1L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)