

Product datasheet for **GA113103**

Human TSC22 domain family, member 4 (TSC22D4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSC22D4 CRISPRa kit - CRISPR gene activation of human TSC22 domain family member 4
Symbol:	TSC22D4
Components:	GA113103G1, TSC22 domain family, member 4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113103G2, TSC22 domain family, member 4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113103G3, TSC22 domain family, member 4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9Y3Q8
Synonyms:	THG-1; THG1; TILZ2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)



Product images:

