

Product datasheet for **GA113104**

Human **STK22C (TSSK3)** activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSSK3 CRISPRa kit - CRISPR gene activation of human testis specific serine kinase 3
Symbol:	TSSK3
UniProt ID:	Q96PN8
Synonyms:	SPOGA3; STK22C; STK22D; TSK3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)