

Product datasheet for **GA113107**

Human TMEM49 (VMP1) activation kit by CRISPRa

Product data:

| | |
|---------------|-------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | VMP1 CRISPRa kit - CRISPR gene activation of human vacuole membrane protein 1 |
| Symbol: | VMP1 |
| UniProt ID: | <u>Q96GC9</u> |
| Synonyms: | EPG3; TANGO5; TMEM49 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)