


Product datasheet for **GA113137**

Human SFXN3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SFXN3 CRISPRa kit - CRISPR gene activation of human sideroflexin 3
Symbol:	SFXN3
Components:	GA113137G1, SFXN3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113137G2, SFXN3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113137G3, SFXN3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BWM7</u>
Synonyms:	BA108L7.2; SFX3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)