

Product datasheet for **GA113145**

Human ISG20L2 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | ISG20L2 CRISPRa kit - CRISPR gene activation of human interferon stimulated exonuclease gene 20 like 2 |
| Symbol: | ISG20L2 |
| UniProt ID: | <u>Q9H9L3</u> |
| Synonyms: | HSD38 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)