


Product datasheet for **GA113147**

Human LAS1L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LAS1L CRISPRa kit - CRISPR gene activation of human LAS1 like, ribosome biogenesis factor
Symbol:	LAS1L
Components:	GA113147G1, LAS1L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113147G2, LAS1L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113147G3, LAS1L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y4W2</u>
Synonyms:	dj475B7.2; Las1-like; WTS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)