

Product datasheet for **GA113148**

Human HYI activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HYI CRISPRa kit - CRISPR gene activation of human hydroxypyruvate isomerase (putative)
Symbol:	HYI
UniProt ID:	Q5T013
Synonyms:	HT036; MGC20767; RP11-506B15.5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)