

Product datasheet for **GA113150**

Human QTRT1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	QTRT1 CRISPRa kit – CRISPR gene activation of human queuine tRNA-ribosyltransferase catalytic subunit 1
Symbol:	QTRT1
Components:	GA113150G1, QTRT1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113150G2, QTRT1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113150G3, QTRT1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BXR0</u>
Synonyms:	FP3235; TGT; TGUT
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

