


Product datasheet for **GA113155**

Human SEH1L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SEH1L CRISPRa kit - CRISPR gene activation of human SEH1 like nucleoporin
Symbol:	SEH1L
Components:	GA113155G1, SEH1L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113155G2, SEH1L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113155G3, SEH1L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96EE3
Synonyms:	SEC13L; Seh1; SEH1A; SEH1B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)