

Product datasheet for **GA113157**

Human ZNF93 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF93 CRISPRa kit - CRISPR gene activation of human zinc finger protein 93
Symbol:	ZNF93
UniProt ID:	P35789
Synonyms:	HPF34; HTF34; TF34; TF34, HPF34, HTF34, ZNF505; zinc finger protein 93; zinc finger protein 505; ZNF505
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)