

Product datasheet for **GA113158**

Human HDHD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HDHD3 CRISPRa kit - CRISPR gene activation of human haloacid dehalogenase like hydrolase domain containing 3
Symbol:	HDHD3
Components:	GA113158G1, HDHD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113158G2, HDHD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113158G3, HDHD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9BSH5
Synonyms:	2810435D12Rik; C9orf158
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)