

Product datasheet for **GA113159**

Human PCDH11Y activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCDH11Y CRISPRa kit - CRISPR gene activation of human protocadherin 11 Y-linked
Symbol:	PCDH11Y
Components:	GA113159G1, PCDH11Y gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113159G2, PCDH11Y gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113159G3, PCDH11Y gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BZA8</u>
Synonyms:	PCDH-PC; PCDH22; PCDHX; PCDHY
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)