

Product datasheet for **GA113163**

Human FCRL4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FCRL4 CRISPRa kit - CRISPR gene activation of human Fc receptor like 4
Symbol:	FCRL4
Components:	GA113163G1, FCRL4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113163G2, FCRL4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113163G3, FCRL4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96PJ5
Synonyms:	CD307d; FCRH4; IGFP2; IRTA1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)