

Product datasheet for **GA113166**

Human SH3BGR13 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SH3BGR13 CRISPRa kit - CRISPR gene activation of human SH3 domain binding glutamate rich protein like 3
Symbol:	SH3BGR13
Components:	GA113166G1, SH3BGR13 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113166G2, SH3BGR13 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113166G3, SH3BGR13 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H299</u>
Synonyms:	HEL-S-297; SH3BP-1; TIP-B1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)