

Product datasheet for **GA113174**

Human **LRRC48 (DRC3)** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	DRC3 CRISPRa kit - CRISPR gene activation of human dynein regulatory complex subunit 3
Symbol:	DRC3
UniProt ID:	<u>Q9H069</u>
Synonyms:	CFAP134; LRRC48
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)