

Product datasheet for **GA113176**

Human **RAB33B** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	RAB33B CRISPRa kit - CRISPR gene activation of human RAB33B, member RAS oncogene family
Symbol:	RAB33B
Components:	GA113176G1, RAB33B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113176G2, RAB33B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113176G3, RAB33B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H082</u>
Synonyms:	SMC2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)