

Product datasheet for **GA113189**

Human PLVAP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLVAP CRISPRa kit - CRISPR gene activation of human plasmalemma vesicle associated protein
Symbol:	PLVAP
Components:	GA113189G1, PLVAP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113189G2, PLVAP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113189G3, PLVAP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BX97</u>
Synonyms:	FELS; gp68; PV-1; PV1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)