

## Product datasheet for **GA113190**

### Human ODAD4 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TTC25 CRISPRa kit - CRISPR gene activation of human tetratricopeptide repeat domain 25
Symbol:	ODAD4
UniProt ID:	<a href="#">Q96NG3</a>
Synonyms:	CILD35
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)