

Product datasheet for **GA113224**

Human FAM167A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FAM167A CRISPRa kit - CRISPR gene activation of human family with sequence similarity 167 member A
Symbol:	FAM167A
Components:	GA113224G1, FAM167A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113224G2, FAM167A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113224G3, FAM167A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96KS9</u>
Synonyms:	C8orf13; D8S265
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)