

Product datasheet for **GA113227**

Human DYNLRB1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DYNLRB1 CRISPRa kit - CRISPR gene activation of human dynein light chain roadblock-type 1
Symbol:	DYNLRB1
Components:	GA113227G1, DYNLRB1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113227G2, DYNLRB1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113227G3, DYNLRB1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NP97</u>
Synonyms:	BITH; BLP; DNCL2A; DNLC2A; ROBLD1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)