

Product datasheet for **GA113230**

Human MS4A8B (MS4A8) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MS4A8 CRISPRa kit - CRISPR gene activation of human membrane spanning 4-domains A8
Symbol:	MS4A8
UniProt ID:	Q9BY19
Synonyms:	4SPAN4; CD20L5; MS4A4; MS4A8B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)