

Product datasheet for **GA113239**

Human **SLC4A9** activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC4A9 CRISPRa kit - CRISPR gene activation of human solute carrier family 4 member 9
Symbol:	SLC4A9
Components:	GA113239G1, SLC4A9 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113239G2, SLC4A9 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113239G3, SLC4A9 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96Q91
Synonyms:	AE4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

